

NOVABONE[®]

BIOACTIVE SYNTHETIC GRAFT

MacroPor Si+

Moldable Continuous Porous Osteostimulative Bone Graft

- Interconnected porosity for accelerated healing
- Formable, robust texture, will not migrate
- Ready to use out of package



NovaBone MacroPor Si+

MacroPor Si+ is as a moldable continuous porous structure that can be fitted to the size and shape of the grafting site. The open structure allows for rapid vascularization and mineralization and permits complete device absorption and replacement by new bone.

MacroPor Si+

| Reference No. | Description |
|---------------|--------------------------|
| NB2501 | 1.6cc, 2.0-5.0mm Morsels |
| NB2502 | 2.6cc, 2.0-5.0mm Morsels |
| NB2505 | 5cc, 2.0-5.0mm Morsels |
| NB2508 | 8cc, 2.0-5.0mm Morsels |
| NB2510 | 10cc, 2.0-5.0mm Morsels |
| NB2515 | 15cc, 2.0-5.0mm Morsels |



Customer Service
855-265-8013

13631 Progress Blvd., #600
Alachua, FL 32615
386-462-7660
Facsimile 386-418-1636