

# NOVABONE<sup>®</sup>

BIOACTIVE SYNTHETIC GRAFT

## NovaBone<sup>®</sup> Particulate


The particulate formulation provides users with a small (90-710 $\mu$ m) granular graft that readily mixes with local bone or blood, forming a slurry that can be packed into a defect.

As new bone generation occurs, the material is resorbed by the body.

ENGINEERED BONE REGENERATION

## Formulated for Small Graft Areas

- Mixes easily
- Exceptional handling and placement
- Ideal for smaller voids and defects



### Novabone Particulate (90 $\mu$ m-710 $\mu$ m)

Reference No.	Description
NB0702	2cc, NovaBone Particulate
NB0705	5cc, NovaBone Particulate
NB0710	10cc, NovaBone Particulate
NB0715	15cc, NovaBone Particulate